



RxAir

Air Purification System



Introducing RxAir!

You may be wondering why you've received a package from UV Flu Technologies that contains information about a device called the ViraTech UV-400. Let us explain.

UV Flu Technologies has been going through some exciting changes over the past several months, and in the process we've made a few updates to our brand. RxAir is the new distribution arm of UV Flu Technologies, Inc. And the RxAir Air Purification System is the same exact device as the UV-400; it just has a different name.

We're more committed than ever to promoting a healthy lifestyle through use of our innovative, patented, FDA-approved ViraTech® air purification technology, and we hope you'll stay tuned to see what else we have in store in the near future.

In the meantime, enjoy breathing in cleaner, healthier air!

Sincerely,
RxAir and UV Flu Technologies

CONTACT INFORMATION

Our customer support team is happy to assist you if you need us, so feel free to reach out!



508-362-5455



info@rxair.com




Customer Support:

RxAir
36 Town Forest Road
Oxford, MA 01540

Corporate Offices:

UV Flu Technologies, Inc.
250 Parkway Drive
Suite 150
Lincolnshire, IL 60069

DESTROYS GREATER THAN 99% OF BACTERIA ON A FIRST-PASS BASIS. CHANGES THE VERY AIR YOU BREATHE.




Here's how it works:

1ST PHASE

The first baffle shows the air flow and creates turbulence

2ND PHASE

Microbes & pollutants killed or inactivated by germicidal ultraviolet light, along with odors and VOCs



KEEPS ON KILLING AND CLEANING

Each time air recirculates, the more contaminants are destroyed... permanently

3RD PHASE

Large particulates are trapped, completing the air purification process

RxAir Unit Specifications:

Power usage:	93 watts
Electrical:	120V/60Hz/1.6A
Unit size:	32"x 12.5" x 9.5"
Weight:	18 lbs.
Safety:	ETL listed
UPC Model :	861642000106

Great for use in:

- Bedrooms
- Nurseries
- Home Offices
- Common Areas
- Health Care Facilities
- Salons
- Restaurants
- Nursing Homes



Using the RxAir

You can begin using your RxAir immediately. Simply plug it in, turn it on and it goes to work. No special wiring needed. Cartridge replacement is quick and easy and is typically only required every 1-4 years depending on use.

Features & Benefits

- ✓ FDA-approved as a Class II Medical Device
- ✓ Features innovative, patented ViraTech® air purification technology
- ✓ Kills greater than 99% of airborne bacteria and reduces concentrations of odors and VOCs
- ✓ Protects those who suffer from respiratory conditions or a weakened immune system
- ✓ Low-maintenance, filterless design
- ✓ Cleans the air up to 3 times per hour in an 800-square-foot room
- ✓ Complements any décor with its sleek, refined design

